











## 2.7 -3.5 -4.0 MM SCREWS

Row	Cat. Code	Barcode	Description	Quantity
1	0301600	8699931010824	UPPER EXTREMITY SCREW TRAY	1
2	20127120027	8698673422742	LOCKING CORTICAL SCREW Ø2,7 MM (12-40mm 2mm incr., 45-60mm 5mm incr.) 2 ea.	38
3	20120120027	8698673413238	CORTICAL SCREW Ø2,7 MM (12-40mm 2mm incr., 45-60mm 5mm incr.) 2 ea.	38
4	20127120035	8698673417571	LOCKING CORTICAL SCREW Ø3,5 MM (12-50mm 2mm incr., 55-60mm 5mm incr.). 12-14-22-24.....60mm 5 ea., 16-18-20mm 10 ea.	125
5	20120120035	8698673417076	CORTICAL SCREW Ø3,5 MM (12-50mm 2mm incr., 55-60mm 5mm incr.). 12-14-22-24.....60mm 3 ea., 16-18-20mm 6 ea.	75
6	20427140040	8698673416826	LOCK.CANCELLOUS SCREW Ø4 MM (14-30mm 2mm incr.) 2 ea.	18
7	20427350040	8698673416918	LOCK.CANCELLOUS SCREW Ø4 MM (35-60mm 5mm incr.) 35-55mm 4 ea. 60mm 2 ea.	22
8	20420140040	8698673422582	CANCELLOUS SCREW Ø4 MM (14-30mm 2mm incr.) 2 ea.	18
9	20420350040	8698673422674	CANCELLOUS SCREW Ø4 MM (35-60mm 5mm incr.) 35-55mm 4 ea. 60mm 2 ea.	22
10	20420160540	8698673422001	CANCELLOUS SCREW Ø4 MM (Semi Thread) (16-30mm 2mm incr.) 2 ea.	16
11	20420350540	8698673422087	CANCELLOUS SCREW Ø4 MM (Semi Thread) (35-60mm 5mm incr.) 4 ea.	24
12	20320735040	8698673429925	MALLEOLAR SCREW Ø4 MM (35-70mm 5mm incr.) 2 ea.	16
13	20820000010	8698673434875	SCREW WASHER Ø10 MM	5



### Recommended Screws:

-  2.7 mm Lock Cortical Screw
-  2.7 mm Cortical Screw
-  3.5 mm Lock Cortical, Screw
-  3.5 mm Cortical
-  4.0 mm Lock Cancellous Screw
-  4.0 mm Cancellous Screw
-  4.0 mm Cancel.Screws (SemiThread)
-  4.0 mm Malleolar Screw
-  10 mm Screw Washer

### LOWPROF. CANCELLOUS SCREW (SEMI THREAD) Ø 4 MM

Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20420160540	8698673422001	20410160540	8698673420304	16	2
20420180540	8698673422018	20410180540	8698673420311	18	2
20420200540	8698673422025	20410200540	8698673420328	20	2
20420220540	8698673422032	20410220540	8698673420335	22	2
20420240540	8698673422049	20410240540	8698673420342	24	2
20420260540	8698673422056	20410260540	8698673420359	26	2
20420280540	8698673422063	20410280540	8698673420366	28	2
20420300540	8698673422070	20410300540	8698673420373	30	2
20420350540	8698673422087	20410350540	8698673420397	35	4
20420400540	8698673422094	20410400540	8698673420403	40	4
20420450540	8698673422100	20410450540	8698673420410	45	4
20420500540	8698673422117	20410500540	8698673420427	50	4
20420550540	8698673477933	20410550540	8698673420434	55	4
20420600540	8698673477940	20410600540	8698673420441	60	4

### MALLEOLAR SCREW Ø 4 MM

Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20320735040	8698673429925	20310735040	8699931015997	35	2
20320740040	8698673429932	20310740040	8699931015980	40	2
20320745040	8698673429949	20310745040	8699931015973	45	2
20320750040	8698673429956	20310750040	8699931015966	50	2
20320755040	8698673429963	20310755040	8699931015959	55	2
20320760040	8698673429970	20310760040	8699931015942	60	2
20320765040	8698673429987	20310765040	8699931015935	65	2
20320770040	8698673429994	20310770040	8699931015928	70	2

### SCREW WASHER Ø 10 MM

Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20820000010	8698673434875	20810000010	8698673453852	10	5

LOCKING CORTICAL SCREW SELFT Ø 2.7 MM					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20127120027	8698673422742	20117120027	8698673449886	12	2
20127140027	8698673422759	20117140027	8698673449893	14	2
20127160027	8698673422766	20117160027	8698673449909	16	2
20127180027	8698673422773	20117180027	8698673449916	18	2
20127200027	8698673422780	20117200027	8698673449923	20	2
20127220027	8698673422797	20117220027	8698673449930	22	2
20127240027	8698673422803	20117240027	8698673449947	24	2
20127260027	8698673422810	20117260027	8698673449954	26	2
20127280027	8698673422827	20117280027	8698673449961	28	2
20127300027	8698673422834	20117300027	8698673449978	30	2
20127320027	8698673422841	20117320027	8698673449985	32	2
20127340027	8698673422858	20117340027	8698673449992	34	2
20127360027	8698673422865	20117360027	8698673450004	36	2
20127380027	8698673437289	20117380027	8698673459045	38	2
20127400027	8698673437296	20117400027	8698673459052	40	2
20127450027	8698673444089	20117450027	8698673459069	45	2
20127500027	8698673444096	20117500027	8698673459076	50	2
20127550027	8698673444102	20117550027	8698673459083	55	2
20127600027	8698673444119	20117600027	8698673459090	60	2

CORTICAL SCREW SELFT Ø 2.7 MM					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20120120027	8698673413238	20110120027	8698673426900	12	2
20120140027	8698673423558	20110140027	8698673426894	14	2
20120160027	8698673423565	20110160027	8698673425958	16	2
20120180027	8698673423572	20110180027	8698673426016	18	2
20120200027	8698673428881	20110200027	8698673426047	20	2
20120220027	8698673428898	20110220027	8698673426078	22	2
20120240027	8698673428904	20110240027	8698673426054	24	2
20120260027	8698673428911	20110260027	8698673425965	26	2
20120280027	8698673428928	20110280027	8698673426023	28	2
20120300027	8698673428935	20110300027	8698673425972	30	2
20120320027	8698673428942	20110320027	8698673426009	32	2
20120340027	8698673428959	20110340027	8698673425996	34	2
20120360027	8698673428966	20110360027	8698673426887	36	2
20120380027	8698673428973	20110380027	8698673425989	38	2
20120400027	8698673428980	20110400027	8698673426870	40	2
20120450027	8698673444027	20110450027	8698673459892	45	2
20120500027	8698673444034	20110500027	8698673459151	50	2
20120550027	8698673444041	20110550027	8698673459168	55	2
20120600027	8698673444058	20110600027	8698673459175	60	2

LOW PROF. CORT. LOCK. SCREW SELFT Ø 3.5 MM					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20127120035	8698673417571	20117120035	8699931033656	12	5
20127140035	8698673417588	20117140035	8699931034035	14	5
20127160035	8698673417595	20117160035	8699931034042	16	10
20127180035	8698673417601	20117180035	8699931034059	18	10
20127200035	8698673417618	20117200035	8699931034066	20	10
20127220035	8698673417625	20117220035	8699931034080	22	5
20127240035	8698673417632	20117240035	8699931034097	24	5
20127260035	8698673417649	20117260035	8699931034103	26	5
20127280035	8698673417656	20117280035	8699931034110	28	5
20127300035	8698673417663	20117300035	8699931034387	30	5
20127320035	8698673417670	20117320035	8699931034394	32	5
20127340035	8698673417687	20117340035	8699931034400	34	5
20127350035	8698673445178	20117360035	8699931034417	36	5
20127360035	8698673417694	20117380035	8699931034424	38	5
20127380035	8698673417700	20117400035	8699931034431	40	5
20127400035	8698673417717	20117420035	8699931034448	42	5
20127420035	8698673417724	20117440035	8699931034455	44	5
20127440035	8698673417731	20117460035	8699931034462	46	5
20127450035	8698673445888	20117480035	8699931034479	48	5
20127460035	8698673417748	20117500035	8699931034561	50	5
20127480035	8698673417755	20117550035	8699931036374	55	5
20127500035	8698673417762	20117600035	8699931034615	60	5

LOW PROF. CORT. SELFT SCREW Ø 3.5 MM					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20120120035	8698673417076	20110120035	8698673419063	12	3
20120140035	8698673417083	20110140035	8698673419070	14	3
20120160035	8698673417090	20110160035	8698673419087	16	6
20120180035	8698673417106	20110180035	8698673419094	18	6
20120200035	8698673417113	20110200035	8698673419100	20	6
20120220035	8698673417120	20110220035	8698673419117	22	3
20120240035	8698673417137	20110240035	8698673419124	24	3
20120260035	8698673417144	20110260035	8698673419131	26	3
20120280035	8698673417151	20110280035	8698673419148	28	3
20120300035	8698673417168	20110300035	8698673419155	30	3
20120320035	8698673417175	20110320035	8698673419162	32	3
20120340035	8698673417182	20110340035	8698673419179	34	3
20120360035	8698673417199	20110360035	8698673419186	36	3
20120380035	8698673417205	20110380035	8698673419193	38	3
20120400035	8698673417212	20110400035	8698673419209	40	3
20120420035	8698673417229	20110420035	8698673419216	42	3
20120440035	8698673417236	20110440035	8698673419223	44	3
20120460035	8698673417243	20110460035	8698673419230	46	3
20120480035	8698673417250	20110480035	8698673419247	48	3
20120500035	8698673417267	20110500035	8698673419254	50	3
20120550035	8698673477216	20110550035	8698673463943	55	3
20120600035	8698673426238	20110600035	8698673419308	60	3

LOW PROF. CANCEL. LOCK. SCREW SELFT Ø 4 MM					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20427140040	8698673416826	20417140040	8698673460546	14	2
20427160040	8698673416833	20417160040	8698673460553	16	2
20427180040	8698673416840	20417180040	8698673460560	18	2
20427200040	8698673416857	20417200040	8698673460577	20	2
20427220040	8698673416864	20417220040	8698673460584	22	2
20427240040	8698673416871	20417240040	8698673460591	24	2
20427260040	8698673416888	20417260040	8698673460607	26	2
20427280040	8698673416895	20417280040	8698673460614	28	2
20427300040	8698673416901	20417300040	8698673460621	30	2
20427350040	8698673416918	20417350040	8698673477612	35	4
20427400040	8698673416925	20417400040	8698673460676	40	4
20427450040	8698673416932	20417450040	8698673460683	45	4
20427500040	8698673416949	20417500040	8698673460690	50	4
20427550040	8698673476639	20417550040	8698673460706	55	4
20427600040	8698673476646	20417600040	8698673460713	60	2

LOWPROF. CANCELLOUS SCREW (FULL THREAD) Ø 4 MM					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20420160040	8698673422599	20410160040	8698673420090	16	2
20420180040	8698673422605	20410180040	8698673420106	18	2
20420200040	8698673422612	20410200040	8698673420113	20	2
20420220040	8698673422629	20410220040	8698673420120	22	2
20420240040	8698673422636	20410240040	8698673420137	24	2
20420260040	8698673426993	20410260040	8698673420144	26	2
20420280040	8698673422650	20410280040	8698673420151	28	2
20420300040	8698673422667	20410300040	8698673420168	30	2
20420350040	8698673422674	20410350040	8698673420175	35	4
20420400040	8698673422681	20410400040	8698673420182	40	4
20420450040	8698673422698	20410450040	8698673420199	45	4
20420500040	8698673422704	20410500040	8698673420205	50	4
20420550040	8698673434394	20410550040	8698673420212	55	4
20420600040	8698673444140	20410600040	8698673420229	60	2