



Lower Extremity

LW-08 Lower Extremity Plates Set

DISTAL FEMUR LATERAL PLATES

Indications: Femoral distal shaft fractures, periprosthetic fractures, supracondylar fractures, intra-articular fractures.

Plate Profile: Thickness: 6.0 mm, Width: 16.0 mm, Distance between holes: 22.4 mm

Cat. Code	Barcode	Description	Side	Material	Length (mm)
35125700109	8698673418240	LOWPRO DIS.FEMUR INT. ANAT.PLATE LAT. 9 HOLES	R	Ti	145
35125700111	8698673418257	LOWPRO DIS.FEMUR INT. ANAT.PLATE LAT. 11 HOLES	R	Ti	188
35125700113	8698673418264	LOWPRO DIS.FEMUR INT. ANAT.PLATE LAT. 13 HOLES	R	Ti	231
35125700115	8698673423268	LOWPRO DIS.FEMUR INT. ANAT.PLATE LAT. 15 HOLES	R	Ti	274
35126700109	8698673418202	LOWPRO DIS.FEMUR INT. ANAT.PLATE LAT. 9 HOLES	L	Ti	145
35126700111	8698673418219	LOWPRO DIS.FEMUR INT. ANAT.PLATE LAT. 11 HOLES	L	Ti	188
35126700113	8698673418226	LOWPRO DIS.FEMUR INT. ANAT.PLATE LAT. 13 HOLES	L	Ti	231
35126700115	8698673427273	LOWPRO DIS.FEMUR INT. ANAT.PLATE LAT. 15 HOLES	L	Ti	274



Recommended Screws:

4.5 mm Cortical Screw

4.5 mm Lock Cortical Screw

6.5 mm Cancellous Screw

6.5 mm Lock Cancellous Screw

