

LW-04 DHS/DCS Plates Set

DCS DYNAMIC CONDYLAR SCREW PLATES



Indications: For supracondylar fractures and intercondylar T and Y type fractures.

Plate Profile: Thickness: 6.5 mm, Width: 19.0 mm, Distance between holes: 16.0 mm

| Cat. Code | Barcode | Description | Material | Length (mm) |
|-------------|---------------|--|----------|-------------|
| 33222950005 | 8698673427754 | LC DYNAMIC CONDYLAR PLATE DCS LOCK 95° 5 HOLES | Ti | 105 |
| 33222950006 | 8698673427761 | LC DYNAMIC CONDYLAR PLATE DCS LOCK 95° 6 HOLES | Ti | 121 |
| 33222950008 | 8698673427785 | LC DYNAMIC CONDYLAR PLATE DCS LOCK 95° 8 HOLES | Ti | 153 |
| 33222950010 | 8698673427808 | LC DYNAMIC CONDYLAR PLATE DCS LOCK 95° 10HOLES | Ti | 185 |
| 33222950012 | 8698673427815 | LC DYNAMIC CONDYLAR PLATE DCS LOCK 95° 12HOLES | Ti | 223 |
| 33222950014 | 8698673427822 | LC DYNAMIC CONDYLAR PLATE DCS LOCK 95° 14HOLES | Ti | 255 |
| 33212950005 | 8698673414167 | DCS PLATE 95° 5 HOLES | S.STEEL | 105 |
| 33212950006 | 8698673414174 | DCS PLATE 95° 6 HOLES | S.STEEL | 121 |
| 33212950008 | 8698673414198 | DCS PLATE 95° 8 HOLES | S.STEEL | 153 |
| 33212950010 | 8698673414068 | DCS PLATE 95° 10 HOLES | S.STEEL | 185 |
| 33212950012 | 8698673414082 | DCS PLATE 95° 12 HOLES | S.STEEL | 223 |
| 33212950014 | 8698673414105 | DCS PLATE 95° 14 HOLES | S.STEEL | 255 |



Recommended Screws

4.5 mm Cortical Screw



4.5 mm Lock Cortical Screw



6.5 mm Cancellous Screw



6.5 mm Lock Cancellous Screw

