

LW-22 Ankle Plates Set


PILON PLATES



Indications: The LCP Pilon Plate is indicated for extra-articular and intra-articular fractures and osteotomies of the distal tibia.

Plate Profile: Thickness: 3.0 mm, Width: 10.0 mm, Distance between holes: 12.0 mm.

Cat. Code	Barcode	Description	Material	Length (mm)
37420700007	8698673463349	LCP PILON PLATE LOCK TI 7 HOLES	Ti	143
37420700009	8698673463356	LCP PILON PLATE LOCK TI 9 HOLES	Ti	167
36720700308	8698673423107	DISTAL TIBIA LOCK. ANTERIOR MALLEOLAR PLATE SHORT	Ti	125
36620700209	8698673423091	DISTAL TIBIA LOCK. MEDIAL MALLEOLAR PLATE SHORT	Ti	143
37410700007	8698673463325	LCP PILON PLATE LOCK STEEL 7 HOLES	S.STEEL	143
37410700009	8698673463332	LCP PILON PLATE LOCK STEEL 9 HOLES	S.STEEL	167
36710700308	8698673449244	DIST.TIBIA ANTER.MALLEOL.LOCK.PLATE STEEL SHORT	S.STEEL	125
36610700209	8698673449091	DIST.TIBIA MEDIAL MALLEOL.LOCK.PLATE STEEL SHORT	S.STEEL	143



PILON PLATE



PILON PLATE
(cut & shaped)



ANTERIOR MALLEOLAR PLATE



MEDIAL MALLEOLAR PLATE

Recommended Screws

2.7 mm Cortical Screw



2.7 mm Lock Cortical Screw



3.5 mm Cortical and Lock Cortical



3.5 mm Lock Cortical, Screw



4.0 mm Cancellous Screw



4.0 mm Lock Cancellous Screw

