

UP-01 Upper Extremity Plates Set

DISTAL HUMERUS EPICONDYLAR LATERAL PLATES



Indications: Intraarticular fractures of the distal humerus, supracondylar fractures and nonunions of the distal humerus.

Plate Profile: Thickness: 3.0 mm, Width: 12.0 mm, Distance between holes: 14.0 mm.

Cat. Code	Barcode	Description	Side	Material	Length (mm)
31825700008	8698673448568	LDC DISTAL HUM. EPICOND.LAT PLATE 8 HOLES	R	Ti	101
31825700010	8698673448575	LDC DISTAL HUM. EPICOND.LAT PLATE 10 HOLES	R	Ti	127
31825700012	8698673448582	LDC DISTAL HUM. EPICOND.LAT PLATE 12 HOLES	R	Ti	157
31825700014	8698673448599	LDC DISTAL HUM. EPICOND.LAT PLATE 14 HOLES	R	Ti	187
31826700008	8698673448612	LDC DISTAL HUM. EPICOND.LAT PLATE 8 HOLES	L	Ti	101
31826700010	8698673448629	LDC DISTAL HUM. EPICOND.LAT PLATE 10 HOLES	L	Ti	127
31826700012	8698673448636	LDC DISTAL HUM. EPICOND.LAT PLATE 12 HOLES	L	Ti	157
31826700014	8698673448643	LDC DISTAL HUM. EPICOND.LAT PLATE 14 HOLES	L	Ti	187
31815700008	8698673448469	LDC DISTAL HUM. EPICOND.LAT PLATE 8 HOLES	R	S.STEEL	101
31815700010	8698673448476	LDC DISTAL HUM. EPICOND.LAT PLATE 10 HOLES	R	S.STEEL	127
31815700012	8698673448483	LDC DISTAL HUM. EPICOND.LAT PLATE 12 HOLES	R	S.STEEL	157
31815700014	8698673448490	LDC DISTAL HUM. EPICOND.LAT PLATE 14 HOLES	R	S.STEEL	187
31816700008	8698673448513	LDC DISTAL HUM. EPICOND.LAT PLATE 8 HOLES	L	S.STEEL	101
31816700010	8698673448520	LDC DISTAL HUM. EPICOND.LAT PLATE 10 HOLES	L	S.STEEL	127
31816700012	8698673448537	LDC DISTAL HUM. EPICOND.LAT PLATE 12 HOLES	L	S.STEEL	157
31816700014	8698673448544	LDC DISTAL HUM. EPICOND.LAT PLATE 14 HOLES	L	S.STEEL	187



Recommended Screws

2.7 mm Cortical Screw



2.7 mm Lock Cortical Screw



3.5 mm Cortical Screw



3.5 mm Lock Cortical Screw



4.0 mm Cancellous Screw



4.0 mm Lock Cancellous Screw

