

LW-08 Lower Extremity Plates Set

PROXIMAL TIBIA MEDIAL PLATES

Indications: The plate is indicated for metaphyseal fractures of the medial tibial plateau and split-type fractures of the medial tibial plateau

Plate Profile: Thickness: 4.0 mm, Width: 15.5mm, Distance between holes: 14.0 mm

| Cat. Code | Barcode | Description | Side | Material | Length (mm) |
|-------------|---------------|---|------|----------|-------------|
| 35025700207 | 8698673418479 | LOWPROF.PROX.TIB.INT.ANAT.PLATE MED. 7 HOLES | R | Ti | 95 |
| 35025700209 | 8698673418486 | LOWPROF.PROX.TIB.INT.ANAT.PLATE MED. 9 HOLES | R | Ti | 123 |
| 35025700211 | 8698673418493 | LOWPROF.PROX.TIB.INT.ANAT.PLATE MED. 11 HOLES | R | Ti | 151 |
| 35025700213 | 8698673418509 | LOWPROF.PROX.TIB.INT.ANAT.PLATE MED. 13 HOLES | R | Ti | 191 |
| 35025700215 | 8698673423480 | LOWPROF.PROX.TIB.INT.ANAT.PLATE MED. 15 HOLES | R | Ti | 231 |
| 35026700207 | 8698673418554 | LOWPROF.PROX.TIB.INT.ANAT.PLATE MED. 7 HOLES | L | Ti | 95 |
| 35026700209 | 8698673418561 | LOWPROF.PROX.TIB.INT.ANAT.PLATE MED. 9 HOLES | L | Ti | 123 |
| 35026700211 | 8698673418578 | LOWPROF.PROX.TIB.INT.ANAT.PLATE MED. 11 HOLES | L | Ti | 151 |
| 35026700213 | 8698673418585 | LOWPROF.PROX.TIB.INT.ANAT.PLATE MED. 13 HOLES | L | Ti | 191 |
| 35026700215 | 8698673440494 | LOWPROF.PROX.TIB.INT.ANAT.PLATE MED. 15 HOLES | L | Ti | 231 |
| 31415000207 | 8698673410602 | ANATOMICPLATE PROX.TIBIA MED. 7 HOLES | R | S.STEEL | 95 |
| 31415000209 | 8698673410619 | ANATOMICPLATE PROX.TIBIA MED. 9 HOLES | R | S.STEEL | 123 |
| 31415000211 | 8698673410589 | ANATOMICPLATE PROX.TIBIA MED. 11 HOLES | R | S.STEEL | 151 |
| 31415000213 | 8698673410596 | ANATOMICPLATE PROX.TIBIA MED. 13 HOLES | R | S.STEEL | 191 |
| 31415000215 | 8698673423497 | ANATOMICPLATE PROX.TIBIA MED. 15 HOLES | R | S.STEEL | 231 |
| 31416000207 | 8698673410640 | ANATOMICPLATE PROX.TIBIA MED. 7 HOLES | L | S.STEEL | 95 |
| 31416000209 | 8698673410657 | ANATOMICPLATE PROX.TIBIA MED. 9 HOLES | L | S.STEEL | 123 |
| 31416000211 | 8698673410626 | ANATOMICPLATE PROX.TIBIA MED. 11 HOLES | L | S.STEEL | 151 |
| 31416000213 | 8698673410633 | ANATOMICPLATE PROX.TIBIA MED. 13 HOLES | L | S.STEEL | 191 |
| 31416000215 | 8698673479548 | ANATOMICPLATE PROX.TIBIA MED. 15 HOLES | L | S.STEEL | 231 |



Recommended Screws:

4.5 mm Cortical Screw

4.5 mm Lock Cortical Screw

6.5 mm Cancellous Screw

6.5 mm Lock Cancellous Screw

