

## ANGLED BLADE PLATES (Child & Infant)

### ANGLED BLADE PLATES (CHILD)

**Indication:** Angled Blade Plate Child is principally used for intertrochanteric femoral osteotomies in children.

**Plate Profile:** These plates have an angle of 90° between blade and shaft. Shaft with 3 dynamic compression holes. Length of blade measured on inner side.

Plate Thickness: 4.0 mm Plate Width: 11.5 mm Hole Spacing: 16.0 mm T-Profile Blade: 4.0 mm - 11.5 mm Material: Stainless Steel

Screws to be used: Fixation with 3.5 mm cortical screws. Shaft with 3 DCP dynamic compression holes.

| Cat. Code   | Barcode       | Description              | Angle | Displacement | Blade Length |
|-------------|---------------|--------------------------|-------|--------------|--------------|
| 32012803503 | 8699931027587 | ANGLED BLADE PLATE CHILD | 90°   | 8 MM         | 35 MM        |
| 32012804003 | 8699931027594 | ANGLED BLADE PLATE CHILD | 90°   | 8 MM         | 40 MM        |
| 32012804503 | 8699931027600 | ANGLED BLADE PLATE CHILD | 90°   | 8 MM         | 45 MM        |
| 32012805003 | 8699931027617 | ANGLED BLADE PLATE CHILD | 90°   | 8 MM         | 50 MM        |
| 32012805503 | 8699931027624 | ANGLED BLADE PLATE CHILD | 90°   | 8 MM         | 55 MM        |

### ANGLED BLADE PLATES (INFANT)

**Indication:** Angled Blade Plate Infant is principally used for intertrochanteric femoral osteotomies in infants up to 5 years of age.

**Plate Profile:** These plates have an angle of 90° between blade and shaft. Shaft with 3 dynamic compression holes. Length of blade measured on inner side.

Plate Thickness: 2.5 mm Plate Width: 10.5 mm Hole Spacing: 10 mm T-Profile Blade: 2.5 mm - 8.5 mm Material: Stainless Steel

Screws to be used: Fixation with 3.5 mm cortical screws. Shaft with 3 DCP dynamic compression holes.

| Cat. Code   | Barcode       | Description               | Angle | Displacement | Blade Length |
|-------------|---------------|---------------------------|-------|--------------|--------------|
| 32012702503 | 8699931027549 | ANGLED BLADE PLATE INFANT | 90°   | 7 MM         | 25 MM        |
| 32012702803 | 8699931027556 | ANGLED BLADE PLATE INFANT | 90°   | 7 MM         | 28 MM        |
| 32012703003 | 8699931027563 | ANGLED BLADE PLATE INFANT | 90°   | 7 MM         | 30 MM        |
| 32012703203 | 8699931027570 | ANGLED BLADE PLATE INFANT | 90°   | 7 MM         | 32 MM        |

