

LC DCP FEMUR PLATES

Indications: Osteosynthesis of bones, such as femur or tibia allowing buttressing, neutralization, tension banding, bridging, compression.

Plate Profile: Thickness: 6.5 mm, Width: 18.5 mm, Distance between holes: 21.5 mm

Cat. Code	Barcode	Description	Material	Length (mm)
30227000006	8698673418714	LC DCP BOW LOCK FEMORAL INT.LOCK. PLATES 6 HOLES	Ti	142
30227000007	8698673418721	LC DCP BOW LOCK FEMORAL INT.LOCK. PLATES 7 HOLES	Ti	168
30227000008	8698673418738	LC DCP BOW LOCK FEMORAL INT.LOCK. PLATES 8 HOLES	Ti	190
30227000009	8698673418745	LC DCP BOW LOCK FEMORAL INT.LOCK. PLATES 9 HOLES	Ti	212
30227000010	8698673418752	LC DCP BOW LOCK FEMORAL INT.LOCK. PLATES 10 HOLES	Ti	233
30227000011	8698673418769	LC DCP BOW LOCK FEMORAL INT.LOCK. PLATES 11 HOLES	Ti	254
30227000012	8698673418776	LC DCP BOW LOCK FEMORAL INT.LOCK. PLATES 12 HOLES	Ti	276
30227000013	8698673418783	LC DCP BOW LOCK FEMORAL INT.LOCK. PLATES 13 HOLES	Ti	298
30227000014	8698673418790	LC DCP BOW LOCK FEMORAL INT.LOCK. PLATES 14 HOLES	Ti	320



Recommended Screws:

4.5 mm Cortical Screw

4.5 mm Lock Cortical Screw

6.5 mm Cancellous Screw

6.5 mm Lock Cancellous Screw