


PROXIMAL HUMERUS LOCKING PLATES

Indications: MISS Proximal Humerus Locking Plate is indicated for fractures and fracture dislocations, osteotomies and nonunions of the proximal humerus, specially for patients with osteopenic bone.

Plate Profile: Thickness: 3.9 mm, Width: 12 mm, Distance between holes: 15 mm

Cat. Code	Barcode	Description	Material	Length (mm)
34620710012	8698673441415	MISS PROXIMAL HUMERUS LOCKING PLATE 12 HOLES	Ti	94
34620710014	8698673441422	MISS PROXIMAL HUMERUS LOCKING PLATE 14 HOLES	Ti	124
34620710016	8698673441439	MISS PROXIMAL HUMERUS LOCKING PLATE 16 HOLES	Ti	150
34620710017	8698673441446	MISS PROXIMAL HUMERUS LOCKING PLATE 17 HOLES	Ti	162
34620710019	8698673469693	MISS PROXIMAL HUMERUS LOCKING PLATE 19 HOLES	Ti	190
34620710021	8698673441453	MISS PROXIMAL HUMERUS LOCKING PLATE 21 HOLES	Ti	218
34610710012	8698673441903	MISS PROXIMAL HUMERUS LOCKING PLATE 12 HOLES	S.STEEL	94
34610710014	8698673441910	MISS PROXIMAL HUMERUS LOCKING PLATE 14 HOLES	S.STEEL	124
34610710016	8698673441927	MISS PROXIMAL HUMERUS LOCKING PLATE 16 HOLES	S.STEEL	150
34610710017	8698673441934	MISS PROXIMAL HUMERUS LOCKING PLATE 17 HOLES	S.STEEL	162
34610710019	8698673478572	MISS PROXIMAL HUMERUS LOCKING PLATE 19 HOLES	S.STEEL	190
34610710021	8698673441941	MISS PROXIMAL HUMERUS LOCKING PLATE 21 HOLES	S.STEEL	218



Recommended Screws

2.7 mm Cortical Screw



2.7 mm Lock Cortical Screw



3.5 mm Cortical Screw



3.5 mm Lock Cortical Screw



4.0 mm Cancellous Screws



4.0 mm Lock Cancellous Screws

