

MISS PROXIMAL TIBIA LATERAL PLATES

Indications: Split-type fractures of the lateral tibial plateau, Lateral split fractures with associated depressions.

Plate Profile: Thickness: 2.0 mm, Width: 10.0 mm, Distance between holes: 12.0 mm.

Cat. Code	Barcode	Description	Side	Material	Length (mm)
35025710107	8698673444836	MISS LC PROXIMAL TIBIA LATERAL PLATE 7 HOLES	R	Ti	103
35025710109	8698673444843	MISS LC PROXIMAL TIBIA LATERAL PLATE 9 HOLES	R	Ti	133
35025710111	8698673438521	MISS LC PROXIMAL TIBIA LATERAL PLATE 11 HOLES	R	Ti	174
35025710113	8698673444850	MISS LC PROXIMAL TIBIA LATERAL PLATE 13 HOLES	R	Ti	215
35026710107	8698673444881	MISS LC PROXIMAL TIBIA LATERAL PLATE 7 HOLES	L	Ti	103
35026710109	8698673444898	MISS LC PROXIMAL TIBIA LATERAL PLATE 9 HOLES	L	Ti	133
35026710111	8698673444911	MISS LC PROXIMAL TIBIA LATERAL PLATE 11 HOLES	L	Ti	174
35026710113	8698673444935	MISS LC PROXIMAL TIBIA LATERAL PLATE 13 HOLES	L	Ti	215
35015700107	8699931023107	MISS LC PROX. TIBIA LATERAL PLATE 7 HOLES	R	S.STEEL	103
35015700109	8699931022544	MISS LC PROX. TIBIA LATERAL PLATE 9 HOLES	R	S.STEEL	133
35015700111	8699931021516	MISS LC PROX. TIBIA LATERAL PLATE 11 HOLES	R	S.STEEL	174
35015700113	8699931021523	MISS LC PROX. TIBIA LATERAL PLATE 13 HOLES	R	S.STEEL	215
35016700107	8699931021530	MISS LC PROX. TIBIA LATERAL PLATE 7 HOLES	L	S.STEEL	103
35016700109	8699931005301	MISS LC PROX. TIBIA LATERAL PLATE 9 HOLES	L	S.STEEL	133
35016700111	8699931024333	MISS LC PROX. TIBIA LATERAL PLATE 11 HOLES	L	S.STEEL	174
35016700113	8699931024340	MISS LC PROX. TIBIA LATERAL PLATE 13 HOLES	L	S.STEEL	215

