

## PROXIMAL FEMUR PLATES

**Indications:** Plate is indicated for the treatment of fractures of the trochanteric region including simple intertrochanteric, reverse intertrochanteric, transverse trochanteric, complex multifragmentary and fractures with medial cortex instability. Proximal femur fractures with ipsilateral shaft fractures. Metastatic proximal femur fractures.

Proximal femur osteotomies, Fractures in osteopenic bone, Nonunions and malunions, Basi/transcervical femoral neck fractures, Subcapital femoral neck fractures, Subtrochanteric femur fractures.

**Plate Profile:** Thickness: 6.5 mm, Width: 18.0 mm, Distance between holes: 21.5 mm

Cat. Code	Barcode	Description	Side	Material	Length (mm)
31625710107	8698673447585	MISS LC PROXIMAL FEMUR PLATE 7 HOLES	R	Ti	155
31625710109	8698673447592	MISS LC PROXIMAL FEMUR PLATE 9 HOLES	R	Ti	198
31625710111	8698673428997	MISS LC PROXIMAL FEMUR PLATE 11 HOLES	R	Ti	241
31625710113	8698673429000	MISS LC PROXIMAL FEMUR PLATE 13 HOLES	R	Ti	285
31626710107	8698673447615	MISS LC PROXIMAL FEMUR PLATE 7 HOLES	L	Ti	155
31626710109	8698673447622	MISS LC PROXIMAL FEMUR PLATE 9 HOLES	L	Ti	198
31626710111	8698673429031	MISS LC PROXIMAL FEMUR PLATE 11 HOLES	L	Ti	241
31626710113	8698673429048	MISS LC PROXIMAL FEMUR PLATE 13 HOLES	L	Ti	285
31615710107	8698673448384	MISS LC PROXIMAL FEMUR LOCK PLATE 7 HOLES	R	S.STEEL	155
31615710109	8698673448391	MISS LC PROXIMAL FEMUR LOCK PLATE 9 HOLES	R	S.STEEL	198
31615710111	8698673448407	MISS LC PROXIMAL FEMUR LOCK PLATE 11 HOLES	R	S.STEEL	241
31615710113	8698673448414	MISS LC PROXIMAL FEMUR LOCK PLATE 13 HOLES	R	S.STEEL	285
31616710107	8698673448421	MISS LC PROXIMAL FEMUR LOCK PLATE 7 HOLES	L	S.STEEL	155
31616710109	8698673448438	MISS LC PROXIMAL FEMUR LOCK PLATE 9 HOLES	L	S.STEEL	198
31616710111	8698673448445	MISS LC PROXIMAL FEMUR LOCK PLATE 11 HOLES	L	S.STEEL	241
31616710113	8698673448452	MISS LC PROXIMAL FEMUR LOCK PLATE 13 HOLES	L	S.STEEL	285



**Recommended Screws:**

4.5 mm Cortical Screw

4.5 mm Lock Cortical Screw

6.5 mm Cancellous Screw

6.5 mm Lock Cancellous Screw

