


SMALL RADIAL HEAD PLATE

Indications: The plate is indicated for extra-articular and intra-articular fractures of the proximal radius and multi-fragmented radial neck fractures.

Plate Profile: Thickness: 1.5 mm, Width: 7.0 mm, Distance between holes: 5.0 mm.

Cat. Code	Barcode	Description	Material	Length (mm)
36920800007	8699931028843	LC SMALL FRAGMENT RADIAL HEAD PLATE 7 HOLES	Ti	33
36920800008	8699931028850	LC SMALL FRAGMENT RADIAL HEAD PLATE 8 HOLES	Ti	38
36920800009	8699931028867	LC SMALL FRAGMENT RADIAL HEAD PLATE 9 HOLES	Ti	43
36910800007	8699931028812	LC SMALL FRAGMENT RADIAL HEAD PLATE 7 HOLES	STEEL	33
36910800008	8699931028829	LC SMALL FRAGMENT RADIAL HEAD PLATE 8 HOLES	STEEL	38
36910800009	8699931028836	LC SMALL FRAGMENT RADIAL HEAD PLATE 9 HOLES	STEEL	43



Recommended Screws:

2.4 mm Lock. Peg Screw

2.7 mm Lock. Peg Screw

2.4 mm Small Head Lock. Screw

2.7 mm Small Head Lock Screw

2.7 mm Small Head Screw



Required Set: UP-03 Small Head Distal Radius Plate Set