

CS-01 HCCS Micro-Mini-Standard Set

HCCS - Micro-Mini-Standard HEADLESS CANNULATED COMPRESSION SCREWS




The innovative design of the screws include a taper, variable thread pitch, and fully threaded length to provide strong compression and holding power for fractures, osteotomies and fusions. These self-drilling screws feature strong guide wires and hex drivers that are designed for optimum screw placement.

Indications;

HCCS Micro: scaphoid fractures, phalangeal fractures, carpal fusions, MCP fusions, radial head fractures, chevron, Akin and Weil osteotomies and OCD repair.
HCCS Mini: scaphoid fractures and nonunions, carpal fusions, radial styloid fractures, MCP fusions, capitellum fractures, bunionectomies: proximal and distal, tarsal fractures, 5th metatarsal fractures, midfoot fusions, OCD repair and osteotomies.
HCCS Standard: talus fractures, malleolar fractures, midfoot fusions, MTP fusions and greater tuberosity fractures.

HEADLESS CANNULATED COMPRESSION SCREWS Micro-Mini-Standard

Cat. Code	Barcode	Description	Diameter (TIP X TAIL)	Length	Qty.
20622528080	8698673465251	MICRO	2.5 X 2.8 MM	8 MM	3
20622528100	8698673465268	MICRO	2.5 X 2.8 MM	10 MM	3
20622528120	8698673465275	MICRO	2.5 X 2.8 MM	12 MM	3
20622528140	8698673465282	MICRO	2.5 X 2.8 MM	14 MM	3
20622528160	8698673465299	MICRO	2.5 X 2.8 MM	16 MM	3
20622528180	8698673465305	MICRO	2.5 X 2.8 MM	18 MM	3
20622528200	8699931028065	MICRO	2.5 X 2.8 MM	20 MM	3
20622528220	8699931028072	MICRO	2.5 X 2.8 MM	22 MM	3
20622528240	8699931028089	MICRO	2.5 X 2.8 MM	24 MM	3
20623338125	8698673424708	MINI	3.3 X 3.8 MM	12.5 MM	3
20623338150	8698673424715	MINI	3.3 X 3.8 MM	15 MM	3
20623338175	8698673424722	MINI	3.3 X 3.8 MM	17.5 MM	3
20623338200	8698673424739	MINI	3.3 X 3.8 MM	20 MM	3
20623338225	8698673424746	MINI	3.3 X 3.8 MM	22.5 MM	3
20623338250	8698673424753	MINI	3.3 X 3.8 MM	25 MM	3
20623338275	8698673424760	MINI	3.3 X 3.8 MM	27.5 MM	3
20623338300	8698673424777	MINI	3.3 X 3.8 MM	30 MM	3
20624045200	8699931003789	STANDARD	4.0 X 4.5 MM	20 MM	3
20624045250	8698673457539	STANDARD	4.0 X 4.5 MM	25 MM	3
20624045300	8698673457546	STANDARD	4.0 X 4.5 MM	30 MM	3
20624045350	8698673457553	STANDARD	4.0 X 4.5 MM	35 MM	3
20624045400	8698673457560	STANDARD	4.0 X 4.5 MM	40 MM	3
20624045450	8699931028881	STANDARD	4.0 X 4.5 MM	45 MM	3
20624045500	8698673457577	STANDARD	4.0 X 4.5 MM	50 MM	3



Row	Cat. Code	Barcode	Description	Quantity
1	00275065020	8698673498006	HCCS-MINI-MICRO-STANDARD SCREW BOX	1
2	02010101004	8699931021868	SOFT SCR. DRIVER CANNULATED QUICK SMALL	1
3	02050100425	8699931005417	SCREW DRIVER CANULLU QUICK Ø 2.5 MM	1
4	02020103025	8699931010930	SCREW DRIVER QUICK TIP 2.5XØ5X130 MM	1
5	02025000170	8698673469396	CANNULATED SCREW DRIVER QUICK HEX.TIP 1.7X80MM	1
6	02000480017	8699931029253	SCREW DRIVER QUICK HEX. TIP 1.7X80 MM	1
7	04551000130	8699931032574	K-WIRE TUBE Ø10XØ8X130 MM	1
8	04551208160	8699931028720	K-WIRE TUBE Ø12XØ8X160 MM	1
9	04551208190	8699931029550	K-WIRE TUBE Ø12XØ8X190 MM	1
10	23410110110	8680858408696	KIRSCHNER WIRE TROCAR POINT 1x110 mm	4
11	23410150110	8698673452909	KIRSCHNER WIRE TROCAR POINT 1X150 MM	4
12	23410160112	8680858408702	KIRSCHNER WIRE TROCAR POINT 1.2x160 mm	4
13	02062011201	8699931028402	CANNULATED DRILL 2X1.1X20X100 MM	1
14	02066251112	8698673496750	CANNULATED DRILL 2.5X1.1X30X120 MM	1
15	02066321212	8698673496767	CANNULATED DRILL 3.2X1.25X50X150mm	1
16	02061900180	8699931008807	LENGTH MEASURING DEVICE - HCCS MIC-MINI-STD.	1
17	03020000115	8698673493353	PENSET 115 MM	1
18	03651000100	8680858408412	K-WIRE SLEEVE FOR MICRO	1
19	03651000200	8680858408429	K-WIRE SLEEVE FOR MINI	1
20	03651000300	8680858408436	K-WIRE SLEEVE FOR STANDARD	1
21	02065300032	8699931008791	HCCS STANDART DRILL SLEEVE Ø 3.2 MM	1
22	02065300025	8699931008784	HCCS MINI DRILL SLEEVE Ø 2.5 MM	1
23	02060001020	8699931031409	HCCS MICRO DRILL SLEEVE 2 MM	1
24	00280136027	8699931029307	DESIGN TRAY 280X136X27 MM	1
	00290145060	8680858424733	CONTAINER 290X145X60 MM	1