

UP-01 Upper Extremity Plates Set

CLAVICLE PLATES






Indications: For all fractures of the clavicle in metaphyseal and diaphyseal areas and for corrective osteotomy.

Plate Profile: Thickness: 3.0 mm, Width: 10.0 mm, Distance between holes: 14.0 mm.

Cat. Code	Barcode	Description	Side	Material	Length (mm)
36625700004	8698673425859	LOW PROF. LOCK. CLAVICLE PLATE 4 HOLES	R	Ti	54
36625700006	8698673425866	LOW PROF. LOCK. CLAVICLE PLATE 6 HOLES	R	Ti	78
36625700008	8698673425873	LOW PROF. LOCK. CLAVICLE PLATE 8 HOLES	R	Ti	98
36626700004	8698673425880	LOW PROF. LOCK. CLAVICLE PLATE 4 HOLES	L	Ti	54
36626700006	8698673425897	LOW PROF. LOCK. CLAVICLE PLATE 6 HOLES	L	Ti	78
36626700008	8698673425903	LOW PROF. LOCK. CLAVICLE PLATE 8 HOLES	L	Ti	98
36615700004	8698673449107	LOCKING CLAVICLE PLATE 4 HOLES	R	S.STEEL	54
36615700006	8698673449114	LOCKING CLAVICLE PLATE 6 HOLES	R	S.STEEL	78
36615700008	8698673449121	LOCKING CLAVICLE PLATE 8 HOLES	R	S.STEEL	98
36616700004	8698673449176	LOCKING CLAVICLE PLATE 4 HOLES	L	S.STEEL	54
36616700006	8698673449183	LOCKING CLAVICLE PLATE 6 HOLES	L	S.STEEL	78
36616700008	8698673449190	LOCKING CLAVICLE PLATE 8 HOLES	L	S.STEEL	98



Recommended Screws

- 2.7 mm Cortical Screw 
- 2.7 mm Lock. Cortical Screw 
- 3.5 mm Cortical Screw 
- 3.5 mm Lock. Cortical Screw 
- 4.0 mm Cancellous Screw 
- 4.0 mm Lock. Cancellous Screw 