


## LW-04 DHS/DCS Plates Set

# DHS DYNAMIC HIP SCREW PLATES

**Indications:** The plate is indicated for intertrochanteric fractures.

**Plate Profile:** Thickness: 6.5 mm, Width: 19.0 mm, Distance between holes: 16.0 mm

Cat. Code	Barcode	Description	Material	Length (mm)
33222350003	8698673427648	LC DYNAMIC HIP SCREW PLATE DHS LOCK 135° 3HOLES	Ti	73
33222350004	8698673427655	LC DYNAMIC HIP SCREW PLATE DHS LOCK 135° 4HOLES	Ti	89
33222350005	8698673427662	LC DYNAMIC HIP SCREW PLATE DHS LOCK 135° 5HOLES	Ti	105
33222350006	8698673427679	LC DYNAMIC HIP SCREW PLATE DHS LOCK 135° 6HOLES	Ti	127
33222350007	8698673427686	LC DYNAMIC HIP SCREW PLATE DHS LOCK 135° 7HOLES	Ti	143
33222350008	8698673427693	LC DYNAMIC HIP SCREW PLATE DHS LOCK 135° 8HOLES	Ti	159
33222350010	8698673427716	LC DYNAMIC HIP SCREW PLATE DHS LOCK 135° 10HOLES	Ti	191
33222350012	8698673427723	LC DYNAMIC HIP SCREW PLATE DHS LOCK 135° 12HOLES	Ti	223
33222350014	8698673427730	LC DYNAMIC HIP SCREW PLATE DHS LOCK 135° 14HOLES	Ti	255
33212350003	8698673414297	DHS PLATE 135° 3 HOLES	S.STEEL	73
33212350004	8698673414303	DHS PLATE 135° 4 HOLES	S.STEEL	89
33212350005	8698673414310	DHS PLATE 135° 5 HOLES	S.STEEL	105
33212350006	8698673414327	DHS PLATE 135° 6 HOLES	S.STEEL	127
33212350007	8698673414334	DHS PLATE 135° 7 HOLES	S.STEEL	143
33212350008	8698673414341	DHS PLATE 135° 8 HOLES	S.STEEL	159
33212350010	8698673414211	DHS PLATE 135° 10 HOLES	S.STEEL	191
33212350012	8698673414235	DHS PLATE 135° 12 HOLES	S.STEEL	223
33212350014	8698673414259	DHS PLATE 135° 14 HOLES	S.STEEL	255



DHS & DCS LAG SCREW					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20720500000	8698673427945	20710500000	8698673414402	50 mm	1
20720550000	8698673427952	20710550000	8698673414419	55 mm	1
20720600000	8698673427525	20710600000	8698673414426	60 mm	1
20720650000	8698673427532	20710650000	8698673414433	65 mm	1
20720700000	8698673427570	20710700000	8698673414440	70 mm	1
20720750000	8698673427587	20710750000	8698673414457	75 mm	1
20720800000	8698673427594	20710800000	8698673414464	80 mm	1
20720850000	8698673427600	20710850000	8698673414471	85 mm	1
20720900000	8698673427617	20710900000	8698673414488	90 mm	1
20720950000	8698673427624	20710950000	8698673414495	95 mm	1
20720000000	8698673427549	20710000000	8698673414365	100mm	1
20720050000	8698673427556	20710050000	8698673414372	105mm	1
20720100000	8698673427563	20710100000	8698673414389	110mm	1

COMPRESSION SCREW					
Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20820250004	8698673427969	20810250004	8698673422322	25 mm	13



**Recommended Screws**

- 4.5 mm Cortical Screw
- 4.5 mm Lock Cortical Screw
- 6.5 mm Cancellous Screw
- 6.5 mm Lock Cancellous Screw

