


TROCHANTER MAJOR STABILIZATION PLATES

Indications: Unstable A2 and A3 fractures, especially multifragmentary ones with a separated greater trochanter. Excessive secondary fracture impact and medialization of the femoral shaft can be avoided by combining the DHS with the trochanter major stabilization plate. Fixation of the greater trochanter with locking head screws reduces varus angulation and limb shortening. Stabilizing the greater trochanter can reduce internal bleeding and restores the function of the gluteus medius. This results in less postoperative pain and a faster mobilization of the patient.

Plate Profile: Thickness: 2.5 mm, Width: 23.5 mm, Distance between holes: 16.0 mm

Cat. Code	Barcode	Description	Material	Length (mm)
3072000008	8699931014211	TROCHANTER MAJOR STABILIZATION PLATE	Ti	132
3071000008	8699931014228	TROCHANTER MAJOR STABILIZATION PLATE	S.STEEL	132
30720000908	8680858409235	TROCHANTER MAJOR STAB.SUPPLEMENTAL PLATE	Ti	35
30710000908	8680858409228	TROCHANTER MAJOR STAB.SUPPLEMENTAL PLATE	S.STEEL	35
30720003908	8680858409259	TROCHANTER MAJOR STAB.SUPPLEMENTAL SMALL PLATE	Ti	25
30710003908	8680858409242	TROCHANTER MAJOR STAB.SUPPLEMENTAL SMALL PLATE	S.STEEL	25
20720000008	8680858407477	SCREW FOR STABILIZATION PLATE	Ti	
20710000008	8680858407484	SCREW FOR STABILIZATION PLATE	S.STEEL	



LOCKING CORTICAL SCREW Ø 3.5 MM

Titanium		Stainless Steel		Length	Qty
Cat. Code	Barcode	Cat. Code	Barcode		
20127160035	8698673417595	20117160035	8699931034042	16	5
20127180035	8698673417601	20117180035	8699931034059	18	5
20127200035	8698673417618	20117200035	8699931034066	20	5
20127220035	8698673417625	20117220035	8699931034080	22	5
20127240035	8698673417632	20117240035	8699931034097	24	5

ADDITIONAL INSTRUMENTS;

Cat. Code	Barcode	Description	Qty
08061008027	8698673495067	LOCKING DRILL GUIDE Ø 2.7X80 MM	2
02057012027	8680858405664	LOCKING DRILL GUIDE Ø 2.7X120 MM	2
00223127200	8699931023114	QUICK DRILL Ø 2.7X200 MM	2
02005000700	8699931028874	DEPTH GAUGE - UPPER EXT. PLATES	1
02020101025	8698673493315	SCREW DRIVER QUICK TIP Ø 2.5XØ5X120 MM	1
02271000001	8699931000115	PIN - PLATE BENDER	2


Recommended Screws

3.5 mm Cortical Screw

4.5 mm Cortical Screw

4.5 mm Lock Cortical Screw

6.5 mm Cannulated Screw



Compatible Instruments Set : LW-04 DHS-DCS Plate Set and CS-06 Cannulated 6.5 mm Screw Set.